


DVC 08-37  
NSN 6910-01-275-4833

**CONVULSANT ANTIDOTE FOR NERVE AGENT  
TRAINING DEVICE (CANATD)**

---



**Functional Description:**

DVC 08-37 functions in the same manner as the CANA autoinjector. It has the same basic appearance as the actual autoinjector, with flanges on the sides and a safety cap to prevent activation. It has a recocking capability and can be operated by individuals in chemical and environmental protective clothing and under conditions of low visibility. Color coding for the training device is blue with white labels and black lettering. After the safety cap is removed and 2-8 pounds of axial pressure is applied, a 3-4 mm plastic tip will protrude, demonstrating activation. The training device does not contain the drug product or needle.

**Purpose of Trainer:**

To aid in the training of individuals for the proper administration of Convulsant Antidote for Nerve Agent (CANA) without the use of the actual injectors.

**Physical Information:**

6.3" L x 1" diameter

**Equipment Required, Not Supplied:**

None

**Special Installation Requirements:**

None

**Power Requirements:**

None

**Training Category/Level Utilized:**

Medical/Level 3

**Training Requirements Supported:**

All MOSCs and Skill Levels

**Applicable Publications:**

None

**DVC 08-37**  
**NSN 6910-01-275-4833**

**Reference Publications:**

FM 8-285

**Source and Method of Obtaining:**

Available through local TSC

**Logistic Responsible Command, Service, or Agency:**

USAMMA